

Yamane technique, the surgical tips and tricks



Rasha Mousa
MD,FRCS
MIOR



The original Yamane technique



What are the tips in pediatric cases ?

- AC maintainer is mandatory
- Tenon ballooning is very easy to happened
- Marfan cases
 - Thin sclera
 - High incidence of RD
- PI



Why do we need to do modifications ?


- To make procedure easier
- Instrumental problems
- Nature of the case
- Surgeon preference
- Avoid complications



The Modified Yamane ,Injector assisted haptic fixation




The Modified Yamane, Using 27 G Trocars



Courtesy of Danielle Cardillo




The Modified Yamane, Using 27 G Trocars




Courtesy of Danielle Cardillo




The Modified Yamane, Leading haptic flanging before extrusion of trailing haptic



The modified Yamane, Temporary suture



The modified Yamane



Take home message

- Yamane technique is a valid option in pediatric patients
- Proper preparation will make it easier
- Any modification that will make the procedure easier for you or decrease complications is always welcomed



Thank you

